



Oxford® MultiFit



Oxford® MultiFit Reflex




Oxford® MultiFit SL



Oxford® MultiFit SL Reflex

## Sizing Guide: Oxford® MultiFit Slings

	Height (cm)	Girth (cm)	Max SWL	Size	Colour Code
<b>Paediatric</b> (XS - S)	75 - 95	45 - 55	227 kg / 35st	1 (XS)	Brown
	95 - 110	50 - 60	227 kg / 35st	2 (XS)	Brown
	110 - 130	55 - 70	227 kg / 35st	3 (S)	Red
	130 - 150	65 - 80	227 kg / 35st	4 (S)	Red
<b>Adult</b> (M - XL)	135 - 160	70 - 85	227 kg / 35st	5 (M)	Yellow
	145 - 170	75 - 90	227 kg / 35st	6 (M)	Yellow
	155 - 180	85 - 100	227 kg / 35st	7 (M)	Yellow
	160 - 190	95 - 115	227 kg / 35st	8 (M)	Yellow
	160 - 190	105 - 130	227 kg / 35st	9 (L)	Green
	160 - 190	125 - 150	227 kg / 35st	10 (L)	Green
	160 - 190	145 - 175	284 kg / 35st	11 (XL)	Blue
	160 - 190	170 - 200	284 kg / 35st	12 (XL)	Blue
<b>MultiFit SL</b> (M - XL)  (See below)	130 - 155	85 - 100	227 kg / 35st	5 (M)	Yellow
	135 - 160	95 - 110	227 kg / 35st	6 (M)	Yellow
	140 - 165	105 - 125	227 kg / 35st	7 (M)	Yellow
	140 - 165	120 - 140	227 kg / 35st	8 (M)	Yellow
	140 - 165	135 - 155	227 kg / 35st	9 (L)	Green
	140 - 165	150 - 175	227 kg / 35st	10 (L)	Green
	140 - 165	170 - 200	284 kg / 35st	11 (XL)	Blue



### Measuring the Client

You will need to take two measurements to find the most suitable sling size for your client; height and girth. See guidance below whilst referring to the diagram and table above for assistance.

#### 1. Height

Take your client's full height as shown. It may prove more accurate to measure height from a recumbent (lying) position.

#### 2. Girth

For female clients, measure around the hips. For male clients and children, measure around the chest just under the arms.

### MultiFit SL

This MultiFit SL incorporates over-sized leg supports to specifically cater for clients with larger than average hips/thighs. Should you find a client's girth size is too large to match the height in the table above, try these for comparison.



**PLEASE NOTE:** The above information is for guidance only. Joerns Healthcare strongly recommends a thorough risk assessment takes prior to lifting/transferring and client to ensure their complete safety, comfort and dignity at all times.