


Sling comfort is everything to do with the correct fit

An ill-fitting sling leads to user discomfort as well as impairing performance.

Our unique anthropometric sizing system is based on the science of universal Vitruvian Man body ratios. Therefore, unlike using generalist sizing templates which can vary between manufacturers, Care & Independence sling sizing is in sync with *individual personal measurements*.

This means that our slings will deliver the correct fit, greatest comfort and optimal performance without the guesswork.



HOW TO MEASURE

Height


For accuracy, measure a person's height whilst they are lying flat on a bed. If you measure nape of neck to coccyx, knee to hip, etc., you will accumulate inaccuracies and almost always result in a person measuring taller than they really are, potentially leading to a sling that is too big for them.

Girth

The second measurement you will need is the girth. This can be done whilst the person is seated or lying down. Measure right around the widest part of the body (usually hips on women and chest just under the arms on men and children).

Sizing Tables

Cross reference height and girth measurements within the appropriate size table to find the best fitting sling size based on your personal measurements.



SLING SIZE GUIDE

Glove Sling® STANDARD	Height cm	Girth cm	SWL* Kg	Size	Binding Colour
Paediatric Unisex	75-95	45-55	230	1	White
	95-110	50-60	230	2	Red
	110-130	55-70	230	3	Yellow
	130-150	65-80	230	4	Orange
Adult Unisex	135-160	70-85	230	5	White
	145-170	75-90	230	6	Red
	155-180	85-100	230	7	Yellow
	160-190	95-115	230	8	Green
	160-190	105-130	230	9	Black
	160-190	125-150	275	10	White
	160-190	145-175	350	11	Red
	160-190	170-200	350	12	Yellow

Glove Sling® FOXY-STYLE	Height cm	Girth cm	SWL* Kg	Size	Binding Colour
Adult Unisex	130-155	85-100	230	5	White
	135-160	95-110	230	6	Red
	140-165	105-125	230	7	Yellow
	140-165	120-140	230	8	Green
	140-165	135-155	230	9	Black
	140-165	150-175	275	10	White
	140-165	170-200	350	11	Red

FOXY STYLE
If a user's girth doesn't co-align with the height column in the standard chart, a Foxy-style sling may offer a better fit.

People with broad hips and larger thighs often find leg supports aren't long enough to comfortably go around their hips and between their thighs. Longer lengths and wider leg sections of the Foxy-design specifically address this issue.

* Safe Working Load. As published in table unless otherwise specified on individual sling data.