

# Medley<sup>®</sup> Ergo Low



Safe for transfer and sleep



---

## Invacare **Medley Ergo Low**

---

The **Medley Ergo Low** from Invacare offers the perfect balance of economy, functionality and quality. The optimised length of the various mattress platform sections provides increased comfort for clients when sitting up in bed. It also reduces shear and friction when the bed sections are being adjusted and prevents clients sliding downwards. The **Medley Ergo Low** offers all the same benefits as the **Medley Ergo**, but boasts extra-low height levels, thus maximising safety for individuals as risk of falling out of bed. The **Medley Ergo Low** can be transformed into the *Invacare Medley Ergo Low Select*. The **Low Select** is a particularly stylish model in natural beech, creating a homely feel suited to any environment.



---

### Optimised dimensions of the mattress platform

---

The mattress platform sections boast optimal dimensions designed in accordance with anthropometric measurements. The **Medley Ergo Low** offers:

- Increased comfort when sitting
- Reduced shear and friction from sliding downwards during profiling
- Suitability for shorter clients

---

### Improved functionality and efficiency

---

The **Medley Ergo Low** has been specially designed for quick transportation to meet the needs of the transient working environments of Community care provision. The **Medley Ergo Low** can be supplied on a transporter kit so the bed can easily be moved and installed. In addition, any assembling can be done without the use of tools or technicians.

---

### Height range supports clients and carers

---

The flexible design of the **Medley Ergo Low** allows the bed to be mounted at a low height, or a higher position if preferred, creating a range of working heights for the carer and transfer heights for the client. This is facilitated by the unique feature of the bed end, which enables the mattress support to be mounted in both upper and lower positions.



**Medley Ergo Low Select**

---

## Features and options

---



---

### Full length side rails

---

Available with wooden or steel full length side rails. Wooden extendable side rails (150mm) are also available.



---

### Scala 2 range

---

A range of collapsible three quarter length steel side rails, available in various heights to cover standard to deep mattress heights.



---

### Safety in and out of bed

---

In the event of a client falling from the bed, a foldable mat can be placed on the floor and attached to the **Medley Ergo Low** for added safety.

**Product code: 1494608**

---

## Options

---



### **Medley Ergo Low**

The height of the **Medley Ergo Low** bed offers increased security for clients enabling them to sleep with real peace of mind. The **Medley Ergo Low** sits low to the floor (210mm in height) which, when combined with a floor mattress, reduces harm from falls. Carers also have the option to raise the bed to a more suitable height when necessary.

The **Medley Ergo Low** makes it easier for clients to get in and out of bed, providing more independence for those who can self transfer.

The **Medley Ergo Low** is dual height adjustable, offering a 210-610mm and 280-680mm height option to suit a wider range of needs.

---

## Features and Options

---



### **Easy to transport**

All **Medley Ergo** beds are easy to assemble and dismantle and can be transported using transporter kits (optional).



### **Lifting Pole**

Can be adjusted in both height and depth (optional).



### **Select design for Medley Ergo**

The homely design of **Medley Ergo Select** can also be retro-fitted (optional).



### **Double height positions for bed ends**

Adjustable height range for **Medley Ergo Low** 210-610mm or 280-680mm.



### **Optimised dimensions of mattress platform**

According to anthropometrical measurements, the optimised dimensions meets the body size of a wider number of clients and reduces shear and friction.



## Features and options



### Mattress support extension

Extension kit for increasing the length of the bed by 150mm.



### Transport fitting

For easy storage and transport.



### Hand controls

Available with or without lockable functions.

## Technical data

For more comprehensive pre-sales information about this product, including the product's user manual, please visit [www.invacare.co.uk](http://www.invacare.co.uk).

	Width outside Width inside	Length outside Length inside	Height adjustment	Mattress support dimensions	Backrest angle	Thigh angle	Legrest lift
<b>Medley Ergo Low</b>	1020 mm 900 mm	2240 mm 2080 mm	280 - 680 mm or 210 - 610 mm	800 - 230 - 300 - 670 mm	0 - 70°	0 - 24°	0 - 11°
	Tilt option <sup>1</sup>	Total weight	Max. user weight	Heaviest part of the product	Complies to the safety standard, IEC 60601-2-52.		
<b>Medley Ergo Low</b>	11°	64 kg	SWL 180 kg 145 kg	17 kg			

1. Please note: the tilt function is only used for clients under medical supervision. If there is a risk of clients using the tilt function without supervision, please use a lockable hand control.

### NB

The bed must not be used by patients under 12 years of age, or by patients with body size equivalent to an average 12 year old or smaller.

## Frame colour and wood work

Light grey

Beech



Please note colours may vary slightly from those displayed above.

## Warranty

3 years frame  
3 years actuator  
2 years handset

The *Invacare Medley Ergo Low* is part of the Safe Patient Handling range.



safepatienthandling

Invacare Limited  
Pencoed Technology Park - Pencoed  
Bridgend CF35 5AQ - United Kingdom  
Tel: +44 1 656 776 222  
Fax: +44 1656 776 220  
E-mail: [uk@invacare.com](mailto:uk@invacare.com)  
Sales Order E-mail: [ordersuk@invacare.com](mailto:ordersuk@invacare.com)  
[www.invacare.co.uk](http://www.invacare.co.uk)

Invacare Ireland Ltd  
Unit 5 - Seatown Business Campus  
Seatown Road - SWORDS - County Dublin - Ireland  
Tel. +353 1 8107084  
Fax +353 1 8107085  
Email: [Ireland@invacare.com](mailto:Ireland@invacare.com)  
[www.invacare.ie](http://www.invacare.ie)

