

Effective Pressure Area Care

MATTRESS OVERLAY

Material:

High grade, antimicrobial TPU material with optimal envelopment and immersion properties to ensure clinical efficacy and patient comfort.

Handwells:

The 'dog bone' handwells facilitate **easy log rolling**, mobilizing and patient transfers.

Indications for use:

- to aid in the prevention of pressure injuries
- can be used for all levels of risk
- to aid in the treatment of all categories of pressure ulcers (I - IV)

Easy Inflation:

- Clear valve for inflation and adjustments.
- Red CPR valve for **rapid deflation**.

Ventilation Holes:

- aiding the even distribution of pressure across the patient's body surface area.

- provides effective ventilation and optimises micro-climate.

- minimises risk of skin maceration as body fluids are drawn away from the patient's skin.

Security Straps:

Corner Straps secure and centre the Mattress Overlay on the bed.

4000A Mattress Overlay



4000 B includes the EHOB Double Action Air Delivery Pump

EHOB® products have been clinically evaluated as part of effective protocol for the prevention and treatment of pressure injuries.

EHOB®