


Technical information



Unistik[®] 3


Designed to provide reassurance

Unistik® 3 puts the sampling experience firmly in your hands, with side-fire activation enabling you to decide exactly when to activate the device and how much pressure to apply. The alignment guide indicates precisely where the device will activate on the patient's finger.



Needle

The hidden needle retracts automatically after use, minimising the risk of re-use, cross-infection and injury.



'Lock-out wings'

Gives a clear visual indication once the device has been discharged, for absolute safety and certainty.



Comfort Zone Technology®



Alignment guide

Indicates precisely where the device will activate.


Side-fire activation

Allows you to decide when to activate the device and how much pressure to apply.



Lot number labelling

Gives you traceability and peace of mind.



Device specifications


Device dimensions	Height	14.0mm
	Width	14.1mm
	Length	51.5mm
Weight	Unistik® 3 Gentle	2.85g
	Unistik® 3 Comfort	2.86g
	Unistik® 3 Normal	2.89g
	Unistik® 3 Extra	2.93g
	Unistik® 3 Neonatal & Laboratory	3.03g
Material	Housing	HIPS plastic
	Lancet	Stainless steel
	Lancet spring	Heat treated stainless steel
Lancet	Unistik® 3 Gentle ●	30G 1.5mm
	Unistik® 3 Comfort ●●	28G 1.8mm
	Unistik® 3 Normal ●●●	23G 1.8mm
	Unistik® 3 Extra ●●●●	21G 2.0mm
	Unistik® 3 Neonatal & Laboratory ●●●●●	18G 1.8mm
	Activation mode	Side fire
	Puncture site	Finger tips (plus heel sampling with Unistik® 3 Neonatal & Laboratory)
	Intended use	Capillary blood sampling

Packing and shipping configuration

The lancets are packaged into boxes, which are then packed into cases for transportation.



Unistik® 3 box



Shipping case

Box quantity	100	200
Boxes per case	40	20
Devices per case	4000	4000
Case weight	13.8kg	13.6kg
Case size	390 x 520 x 290mm	297 x 587 x 315mm

Ordering information

Contact details:

David Ramsey +44 (0) 1993 812021

United Kingdom, Belgium, Bulgaria, Croatia, Czech Republic, Denmark, Estonia, Finland, Hungary, Ireland, Latvia, Lithuania, Luxembourg, Poland, Netherlands, Romania, Slovakia, Slovenia, Iceland, Norway, Sweden and Liechtenstein.

Thomas Rossbach +49 (0) 6026 97750

Germany, Austria and Switzerland.

Eric Del Cotto + (33) 2 32 51 88 70

France, Italy, Greece, Malta, Portugal and Spain.

Ordering lead times:

Please contact your local Owen Mumford representative for information.

Pricing:

Please contact your local Owen Mumford representative for information.

	Box quantity*	Order code
Unistik® 3 Gentle	100	AT1022
	200	AT1024
Unistik® 3 Comfort	100	AT1042
	200	AT1044
Unistik® 3 Normal	100	AT1002
	200	AT1004
Unistik® 3 Extra	100	AT1012
	200	AT1014
Unistik® 3 Neonatal & Laboratory	100	AT1052
	200	AT1054

* Other box quantities are also available. Please contact your local Owen Mumford representative for information.