

Efficacy details

| Effective against | Test | Contact / Kill time | |
|--|----------|---------------------|------------------|
| | | Dirty conditions | Clean conditions |
| Bactericidal | | | |
| Enterococcus hirae | EN 16615 | 1 minute | 1 minute |
| | EN 13727 | 1 minute | 30 seconds |
| Pseudomonas aeruginosa | EN 16615 | 1 minute | 1 minute |
| | EN 13727 | 1 minute | 30 seconds |
| Staphylococcus aureus | EN 16615 | 1 minute | 1 minute |
| | EN 13727 | 1 minute | 30 seconds |
| Escherichia coli | EN 13727 | 1 minute | 30 seconds |
| | EN 13727 | 1 minute | 30 seconds |
| Mycobactericidal | | | |
| Mycobacterium avium | EN 14348 | | 1 minute |
| Mycobacterium terrae | EN 14348 | | 1 minute |
| Virucidal | | | |
| Adenovirus type 5 | EN 14476 | 2 minutes | 1 minute |
| | EN 16777 | | 1 minute |
| Murine Norovirus | EN 14476 | 2 minutes | 1 minute |
| | EN 16777 | | 1 minute |
| Poliovirus | EN 14476 | | 2 minutes |
| Vaccinia virus | EN 14476 | 1 minute | 1 minute |
| Coronavirus* | EN 14476 | 30 seconds | 1 minute |
| *Tested against Human coronavirus, strain 229E | EN 16777 | | 1 minute |
| Yeasticidal and Fungicidal | | | |
| Aspergillus brasiliensis | EN 13624 | | 2 minutes |
| Candida albicans | EN 16615 | 1 minute | 1 minute |
| | EN 13624 | 1 minute | 1 minute |

Additional Organisms

| | | | |
|---|----------|----------|------------|
| Bactericidal | | | |
| Enterococcus faecium | EN 13727 | | 30 seconds |
| Acinetobacter baumannii | EN 13727 | 1 minute | 30 seconds |
| Klebsiella pneumoniae | EN 13727 | 1 minute | 30 seconds |
| MRSA (Methicillin Resistant Staphylococcus aureus) | EN 13727 | | 1 minute |
| VRE (Vancomycin Resistant Enterococcus faecalis) | EN 13727 | | 1 minute |
| MRSE (Methicillin Resistant Staphylococcus epidermidis) | EN 13727 | | 2 minutes |
| Multidrug Resistant Acinetobacter baumannii | EN 13727 | | 2 minutes |
| VRE (Vancomycin Resistant Enterococcus faecium) | EN 13727 | | 2 minutes |
| Carbapenem Resistant Escherichia coli | EN 13727 | | 2 minutes |
| Carbapenem Resistant Klebsiella pneumoniae | EN 13727 | | 2 minutes |
| Multidrug Resistant Pseudomonas aeruginosa | EN 13727 | | 2 minutes |
| Yeasticidal | | | |
| Candida Auris | EN 13624 | | 1 minute |

Pal Disinfectant Wipes have been tested to EN standards under clean and dirty conditions using the wipe or liquid squeezed from the wipe.

Disinfectant V15 Feb 2024