FEATURED PRODUCT









CODE PRODUCT706547 Integral Mobile Air Care Chair

MOBILE AIR CHAIR

PRESSURE RELIEF SEATING SOLUTION

Adjustable Foam Matrix

A matrix of air flow chambers promotes skin micro-climate and evenly redistributes pressure. A higher-density outer perimeter improves midline sitting, while foam cores are movable for customised positioning.

Patient Surface

Contoured air bladder in the seating surface can be easily inflated and adjusted with the integrated hand pump.

Side Transfer Compatible

Allows for easy side patient transfers for improved WH&S.

Ergonomic Attendant Grip Handles

Control levers for tilt & recline, push handle for attendant use.

Lockable Castors with Directional Control

Front directional castor improves the ability to manoeuvre in a straight line through corridors, making it safer and lighter for attendant use.

Specification:

- 3 year Frame Warranty
- 1 Year Fabric Warranty
- 180kg Safe Working Load
- Choice of 410 or 500mm Seat Widths
- Comes with Padded Footrests
- Genuine Dartex[™] PU Fabric which is anti-bacterial/ anti-microbial and waterproof yet breathable.

TURN OVER FOR MORE DETAILS...



FEATURED PRODUCTS 2019

PRESSURE RELIEF SEATING SOLUTION





TRUE LAY FLAT POSITIONING

Unique patented design ensures true lay flat positioning whilst SmartTilt intelligently manages chair positioning to prevent unintentional head down Trendelenburg.

CLINICAL CONSIDERATIONS

Pressure Injury Management

 Pressure injury management and prevention is the cornerstone of the Mobile Air Chair's design

Optimal Positioning and Posture

- Gentle contouring in the sacral region assists the user in achieving and maintaining optimal positioning and a midline posture
- In the backrest, a higher-density foam outer perimeter improves midline sitting, while foam cores are movable for customised positioning needs
- Swing-away backrest wings and armrests decrease the demands on carers when transferring and positioning the user

Fatigue Management

 Integrated tilt, recline and leg raise allow the Mobile Air Chair to achieve bed-like positioning to facilitate fatigue management and relaxed positioning

Oedema Management

 The generous 90° leg raise angle assists with the management of lower limb oedema

Safe Handling

 Ergonomic grip handles assist with leverage, while a directional locking castor makes manoeuvring in a straight line through long corridors easy and safe

Contact 01392 823233 or email mail@clhgroup.co.uk to request a free no-obligation demonstration in your home.





