

## FEATURED PRODUCTS



## TREAT-EEZI

A technological breakthrough in pressure care, helping to heal stages 2, 3 & 4 bed/pressure sores.

**'Deceptively simple yet incredibly effective'.** **Not to be compared with anything else on the market today,** this patented and NHS clinically approved product is rapidly becoming the preferred purchase by private individuals throughout the community. This is the softest, most easy to use and cost-effective way to treat patients susceptible to pressure sores, allowing the patient a peaceful night's sleep **ANYWHERE!**

Made from thousands of polyester fibres in multiple layers, the Treat-Eezi surface conforms to the natural body contours in both supine and seated positions. This provides pressure relief well below that of normal capillary blood flow, hence skin breakdown and shearing plus friction is simply diminished.

### Features of the Treat-Eezi include:

- The latest and easiest solution to pressure sores
- Keeps pressure below 28 millimetres of mercury (mmHg) on any given part of the body
- Wicks away moisture - reducing the problems caused by sweating
- Breathable - helps regulate body temperature
- Easy to keep clean - machine washable to 80°C
- Can be washed and bleached
- Totally eliminates shearing forces
- Minimises friction between skin and pad
- Easily transported
- No pump - maintenance free and not dependant on electricity
- 2 years warranty



VIEW OVERLEAF FOR PRICING...



For expert information and advice , call us today on 01392 823233  
Alternatively, visit [www.clhgroup.co.uk](http://www.clhgroup.co.uk) or email [info@clhgroup.co.uk](mailto:info@clhgroup.co.uk)



### 1. COMMUNITY PRESSURE ULCER OVERLAY

Available in 3 or 4 layer of fibres for different weight capacities, the Community bed pad will cover part of the bed and will protect the most vulnerable areas. It is ideally suited for switching to the chair when the patient is sitting during the day.



### 2. FULL LENGTH OVERLAY

4 layers of fibres for maximum pressure relief covering the full mattress.



### 3. PILLOWCASE

2-layer surface for pressure relief on and around the head area and especially around the ears.



### 4. SEAT PAD

For use in conjunction with a standard 2, 3 or 4" foam cushion, this 4-layer surface transforms the seat into a high-risk product.



### 5. VAPOUR PERMEABLE COVERS

Ideal for incontinent users and for complying with infection control when switching users. Flaps cover zips and crib 5 tested.

Code	Description	Size	Weight Capacity
706082	Community Pad 3 Layer	100 x 80cm	17 stone
706083	Community Pad 4 Layer	100 x 80cm	25 stone
706081	Full Length Overlay 4 Layer	195 x 90cm	25 stone
706084	Pillowcase 2 Layer	50 x 70cm	11 stone
702593	Seat Pad 4 Layer	50 x 40cm	25 stone
706086	Community Pad Cover	100 x 80cm	n/a
706085	Full Length Overlay Cover	195 x 90cm	n/a
702594	Seat Pad Cover	50 x 40cm	n/a

**Call us today on 01392 823233 for expert information and advice**

Prices exclude VAT

