# GP DETERGENT CLOVE



## CONCENTRATED LIQUID DETERGENT

- Mild ideal for hand use
- Leaves washed items sparkling and grease free
- Long lasting foam

### USES

A detergent concentrate containing additives for efficient grease removal.

A thick, green liquid consisting of detergents, grease-dissolving surfactants, foam boosters, preservative and dye. GP detergent is pH-balanced for mildness to hands.

Designed for cleaning glassware, cutlery, crockery, utensils and worktops, GP Detergent is ideal for use on aluminium, glass, stainless steel, vinyl, plastic, ceramic, wood and enamel.

Code: 418




#### DIRECTIONS

#### DISHWASHING

Dilute 1 part GP Detergent with between 750 to 1000 parts warm water in a washing-up bowl. Add items to be washed and allow a few minutes to soak. Clean items using a sponge, brush or pad.


#### **GENERAL CLEANING**

Dilute 1 part GP Detergent with 300 parts cold or warm water.



#### **DILUTIONS**

GP DETERGENT : WATER	
Dishwashing	General Cleaning
1:750 - 1:1000	1:300



5L