

ALERTA FALL PREVENTION BROCHURE

CLH have been supplying fall prevention alarm mats into care homes for over 10 years and have built up vast knowledge and experience in this sector. We have two main ranges of Alerta fall prevention products; the all new **wireless** and the traditional **wired**. This guide will enable you to easily pick out the products that suit your requirements.

Don't hesitate to contact us for expert advice on how these products can integrate with your existing alarm system. With many different adaptors we are able to provide a solution for most systems.

ALERTA WIRELESS RANGE



The Alerta Wireless Range removes the need for trailing cables from your fall prevention equipment to your nurse call point.

- Enables floor, bed and chair alarm mats as well as nurse call buttons and motion sensors to be connected wirelessly to the nurse call system.
- Can also function as a stand-alone system, giving localised alerts.
- Uses radio frequency technology.
- Range of up to 95 meters.
- Once paired with Alerta devices, it will only transmit to those devices and will not interfere with other Alerta Wall Point Receivers or devices in use nearby.

ALERTA WIRED RANGE



The Alerta Wired Range is our standard range of fall prevention equipment.

- Designed to be plugged into the nurse call point.
- Chair and bed devices come with their own alarm monitor enabling them to function as a standalone unit giving localised alerts.
- The monitor can be powered by mains (adaptor sold separately) or battery making the system extremely versatile.
- Many different adaptors available to ensure these mats work correctly with your system.

For expert advice, call 01392 823233 or email info@clhgroup.co.uk



WIRELESS WALL POINT RECEIVER

The Alerta Wireless Wall Point Receiver connects into the nurse call point on the wall, through a short cable connection, allowing wireless communication with all wireless devices.

- Removes trip hazards of trailing wires.
- Greatly reduces breakages or damage to the cables as on wired mats.
- Allows connection of up to 5 devices to 1 receiver.
- Battery powered with low battery alert.
- For nursing purposes, it can be muted for 3 minutes before automatically resetting.

Code	Description	CLH Price
MULTIPLE	Alerta Wireless Wall Point Receiver	£90.00

WIRELESS DETECT MOTION SENSOR



The Alerta Detect uses infrared technology to trigger an alert when body movement is detected within the range of the infrared beam. It is most effective when mounted at the door or bedside.

Connects wirelessly to the wall point receiver or can be wired directly into the nurse call system.



WIRELESS NURSE CALL BUTTON



The Wireless Nurse Call Button provides the user with assurance that help is at their fingertips. It comes with a lanyard and a holder that can be wall mounted.

Connects wirelessly to the wall point receiver it has been paired with.

Code	Description	CLH Price	Code	Description	CLH Pric
701210	Alerta Wireless Detect Motion Sensor	£77.50	701209	Alerta Wireless Nurse Call Button	£35.00
701164	Alerta Detect Sensor - 2 Pin Mono	£77.50			
701165	Alerta Detect Sensor - 3 Pin Mono	£77.50			

View the full range online at www.clhgroup.co.uk

WIRELESS CHAIR ALERTAMAT

The wireless chair mat senses that the occupant has left the chair and transmits a signal via its inbuilt transmitter to the Wall Point Receiver (effective up to 95 metres away). This provides the caregiver with advanced warning when a user leaves their chair. When the user returns to the chair the alarm is reset automatically.

Comes with transmitter that connects wirelessly to the wall point receiver it has been paired with.

Code	Description	CLH Price
701197	Wireless Chair Alertamat	£90.00

WIRELESS BED ALERTAMAT

The wireless bed mat senses that the occupant has left the bed and transmits a signal via its inbuilt transmitter to the Wall Point Receiver (effective up to 95 metres away) to provide the caregiver with advanced warning when a user at risk of falling leaves their bed. When the user returns to the bed the alarm is reset automatically.

Comes with transmitter that connects wirelessly to the wall point receiver it has been paired with.

Code	Description	CLH Price
701206	Wireless Bed Alertamat	£100.00

E100.00

alert

Dimensions: 25 x 76cm (10 x 30")

Dimensions: 25 x 38cm (10 x 15")

But Milida

E90.00

WIRELESS FLOOR ALERTAMAT

The wireless floor mat is most effective when placed at the bedside or doorway. When it senses pressure on its surface it transmits a signal via its inbuilt transmitter to the Wall Point Receiver (effective up to 95 metres away) to provide the caregiver with advanced warning when a patient who is at risk has left their safe zone.

Comes with transmitter that connects wirelessly to the wall point receiver it has been paired with.

Code	Description	CLH Price
701207	Wireless Floor+ Alertamat	£103.50
701208	Wireless Deluxe Floor Alertamat	£135.00



Dimensions: 60 x 90cm (23.5 x 35.5")

For further information, call 01392 823233 or email info@clhgroup.co.uk



WIRED RANGE

Adaptor cables available to suit the following nurse call systems: Courtney Thorne, Aidcall, Medicare, Quantec, SAS Network II, Eclipse, Touchsafe Pro, Intercall and many more.

CHAIR ALARM MAT & MONITOR

This mat comes with a monitor that enables it to be used as a standalone unit, or through a nurse call system using the lead provided. The mat alarms when the patient leaves the chair.

Code	Description	CLH Price
701118	Alerta Chair Mat - 2 Pin Mono	£95.00
701130	Alerta Chair Mat - 3 Pin Stereo	£95.00

BED ALARM MAT & MONITOR

This mat comes with a monitor that enables it to be used as a standalone unit, or through a nurse call system using the lead provided. The mat alarms when the patient leaves the bed.

Code	Description	CLH Price
701109	Alerta Bed Mat - 2 Pin Mono	£95.00
701121	Alerta Bed Mat - 3 Pin Stereo	£95.00

FLOOR ALARM MAT

These mats wire straight into the nurse call system and alarm when they are stepped on.

Code	Description	CLH Price
701030	Alerta Floor+ Mat - 2 Pin Mono	£69.50
701031	Alerta Floor+ Mat - 3 Pin Stereo	£69.50
701185	Alerta Wood Floor Mat - 2 Pin Mono	£99.50
701186	Alerta Wood Floor Mat - 3 Pin Stereo	£99.50
701101	Deluxe Carpet Mat - 2 Pin Mono	£105.00
701102	Deluxe Carpet Mat - 3 Pin Stereo	£105.00

NURSE CALL BUTTON

This button wires directly into the call alarm system and triggers the alarm when pressed.

Code	Description	CLH Price
701173	Alerta Call Button - 2 Pin Mono	£16.50
701174	Alerta Call Button - 3 Pin Stereo	£16.50











All prices shown exclude VAT



CLH Group Ltd. Registration Number: 10026607 (England & Wales) VAT Number: 236 3082 23 Registered Office: Devonshire House, Marsh Barton Trading Estate, Exeter, Devon EX2 8QW